

# SPECIFICATIONS

CapSureFix Novus MRI™ SureScan™  
5076 Lead

## MR Conditional



Initial Implant  
Existing Implant

## Physical Characteristics

Polarity	Bipolar
Shape	Straight
Chamber	Ventricle or Atrium
MR-Conditional Lengths	35, 45, 52, 58, 65, 85 cm
Connector	IS-1 BI

## Materials

Insulator	Silicone (MED-4719) (inner and outer)
Conductor	MP35N Nickel Alloy
Helix Electrode	Platinum Alloy with Titanium Nitride coating
Ring Electrode	Platinum Alloy with Titanium Nitride coating
Connector Ring	Stainless steel
Connector Pin	Stainless steel

## Diameter

Body	2.0 mm (6 Fr)
------	---------------

## Recommended Introducer Size

Without Guide Wire	7 Fr
With Guide Wire	9 Fr

## Electrode Configuration

Extendable/Retractable Helix Screw
------------------------------------

## Electrode Surface Area

Helix	4.2 mm <sup>2</sup>
Ring	22 mm <sup>2</sup>



## Other Characteristics

Steroid	Dexamethasone acetate
Tip-to-Ring Spacing	10 mm
Unipolar Conductor Resistance	33 $\Omega$ (58 cm)
Bipolar Conductor Resistance	61 $\Omega$ (58 cm)
Serial Prefix	PJN

## Stylets

Inserted	Gray Straight
Packaged	Gray Straight Blue Straight Gray J-shaped* Blue J-shaped* White J-shaped*

\*J-shaped stylets are not available for leads 65 or 85 cm.

## Brief Statement

See the device manual for detailed information regarding the instructions for use, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. If using an MRI SureScan® device, see the MRI SureScan® technical manual before performing an MRI. For further information, contact your local Medtronic representative and/or consult the Medtronic website at [medtronic.eu](http://medtronic.eu).

For applicable products, consult instructions for use on [www.medtronic.com/manuals](http://www.medtronic.com/manuals). Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

**Important Reminder:** This information is intended only for users in markets where Medtronic products and therapies are approved or available for use as indicated within the respective product manuals. Content on specific Medtronic products and therapies is not intended for users in markets that do not have authorization for use.

# Medtronic

## Europe

Medtronic International Trading Sàrl.  
Route du Molliau 31  
Case postale  
CH-1131 Tolochenaz  
[www.medtronic.eu](http://www.medtronic.eu)  
Tel: +41 (0)21 802 70 00  
Fax: +41 (0)21 802 79 00

## United Kingdom/Ireland

Medtronic Limited  
Building 9  
Croxley Park  
Hatters Lane  
Watford  
Herts WD18 8WW  
[www.medtronic.co.uk](http://www.medtronic.co.uk)  
Tel: +44 (0)1923 212213  
Fax: +44 (0)1923 241004

UC201305283aEE © Medtronic 2020.  
All rights reserved.